


Description	
Source	HEK293 Met1-Pro281, with a C-terminal DDK and His tag Accession # NP_112230
Purity	>90% by SDS-PAGE
N-terminal Sequence	Arg26 predicted
Predicted Molecular Mass	29.2kDa
SDS-PAGE	40kDa, reducing conditions
Formulation	Lyophilized from 1X PBS with Trehalose
Preparation and Storage	
Reconstitution	Reconstitute at 50 µg/mL in sterile PBS containing at least 0.1% protein
Shipping	Ambient or cold
Storage	5 years at -80°C; 2 years at -20°C or 4°C Upon reconstitution: 5 years at -80°C; 2 years at -20°C
Expiration	1/28/2021



CTRP1 on a Coomassie stained 4-12% SDS-PAGE gel, reducing conditions.