


Description	
Source	HEK293 Met1-Lys414, with a C-terminal DDK and His tag Accession # NP_776249
Purity	>95% by SDS-PAGE
N-terminal Sequence	Leu21, predicted
Predicted Molecular Mass	45.1kDa
SDS-PAGE	56kDa, reducing conditions
Formulation	Lyophilized from 1X PBS with Trehalose
Preparation and Storage	
Reconstitution	Reconstitute at 50 µg/mL in sterile PBS containing at least 0.1% protein
Shipping	Ambient
Storage	5 years at -80°C; 2 years at -20°C or 4°C Upon reconstitution: 5 years at -80°C; 2 years at -20°C
Expiration	06/2022



Vaspin on a Coomassie stained 4-12% SDS-PAGE gel, reducing conditions.